

Hill-Rom® 900 Accella™ bed

Improve safety. Enhance outcomes. Simplify tasks.



Features and benefits

Improve patient safety

Reduce patient discomfort with **SlideGuard™**, especially designed to minimise the need for patient repositioning

Caregivers can respond more rapidly with **Auto CPR**, which flattens the back rest and automatically places the bed into a horizontal position, should patient complications occur

Alerts caregivers of potential bed exit via the **3-Mode Bed Exit Alarm**, available with remote monitoring capabilities*

Enhance patient outcomes

One-Touch Side Egress positions the bed for safer and easier bed exit

Stable support is provided during bed exit via an **ergonomically-placed side egress grip handle**

Optimise patient care with **Head-of-Bed Angle Alarm**, alerting caregivers of any safety concerns

Optional **Accella Therapy™ integrated surface** for enhanced protection against skin breakdown

* Available in selected markets





Accella Therapy integrated surface controls



3-Mode Bed Exit Alarm



Auto CPR



Single connection in the bed frame to easily connect the Accella Therapy surface to the bed



One-Touch Side Egress positioning

Features and benefits

Simplify caregiver workflow

Simple Weigh improves workflow efficiency with an integrated scale and BMI calculation, excluding the weight of hanging accessories and IV poles, therefore reducing the steps needed to obtain patient weight

Available on the right or left side of the bed, the bright and easy-to-use **Graphical Caregiver Interface (GCI)** offers straightforward visuals of patient information, such as head-of-bed angle and other caregiver alarms

Intuitive bed controls consistent with other Hill-Rom beds, requiring minimal in-service staff training helping to mitigate the risk of user error

Optional in-bed x-ray, via a selection of Hill-Rom therapeutic surfaces and frame accessory



Simple Weigh in-bed scale

Technical Specifications		
Safe Working Load		250 kg
Bed weight		170 kg
Height adjustment		38.6 cm to 80.8 cm
SlideGuard sliding backrest		up to 12 cm
Protection class		IPX4
Bed angles	Head section	0 to 65°
	Thigh section	0 to 28°
	Foot section	-22° to -3°
	Trendelenburg/Reverse Trendelenburg	-17 to +17°
Sleep deck	Width	90 cm
	Length	200 cm
Siderail height		39.3 cm
Simple Weigh (Class III OIML)		●
Body Mass Index calculation		●
Bed extension with linen holder		●
Audible brake-off indicator		●
Intelligent night light		●
Accella Therapy connection ready solution		●
Controls	Integrated in siderails	●
	Pendant	○
	Bilateral high-low foot control	○
Graphical Caregiver Interface (one side of the bed only)	Integrated in right siderail	●
	Integrated in left siderail	○
Steering castor	At head end	●
	At foot end	○
Castors	150 mm integral	●
	125 mm double band	○
	150 mm double band	○
5th castor		○
Brake/steer pedal at head and foot		○
3-Mode Bed Exit Alarm		○
Wired connection for 3-Mode Bed Exit Alarm		○
Nurse call		○

Additional standard features for all configurations

- 6 attachment points for patient restraints
- Lifting/IV pole holders
- Unweighed bilateral drainage hanging bar with hooks
- Electric cord holder
- Adjustable mattress retainers
- Battery backup
- Bilateral Mechanical CPR lever
- Removable head/footboards
- Head-of-Bed Angle Alarm
- EasyChair and flat position function
- One-Touch Side Egress
- Boost function
- High Pressure Laminate (HPL) sleep deck surface
- Central brake at foot end
- Low Height Indicator
- Electric CPR
- Service required indicator

● = standard ○ = optional

Hill-Rom is a leading global medical technology company whose products, services and more than 10,000 employees worldwide help people get better care inside and outside the hospital. Our innovations in five core areas – **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency, and Respiratory Health** – improve clinical and economic outcomes and ensure caregivers in more than 100 countries have the products they need to protect patients, speed up recoveries and manage conditions. **Every day, around the world, we enhance outcomes for patients and their caregivers.** Learn more at hill-rom.com.

mobility is life 

Hill-Rom is committed to the concept that mobilising people early and supporting their independence, improves life - in the hospital and at home.

The Hill Rom 900 Accella bed is intended for use in intensive, acute and ambulatory care settings for patients with a weight ≥ 40 kg, height ≥ 146 cm and BMI ≥ 17 .

Medical device (93/42/EEC): Class I

Weigh system (2014/31/EEC): Class III

Manufacturer: Hill-Rom S.A.S. – BP 14 ZI du Talhouët 56330 Pluvigner France

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2018 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. Doc. No: 5EN134301-03, 06 June 2018

Not all products/options are available in all countries.

For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

Enhancing outcomes for
patients and their caregivers:

Hill-Rom